


How to pair signature pad with Concentric - It's a 2 step process


Step 1 - Configure

- a** Plug the signature pad into the workstation, log in to Concentric and click on icon on right hand side (your initials).



- b** Select 'Configure signature pads'.

- c** Then the below message will appear, click on 'Pair a signature pad'



Configure signature pad

i A Wacom STU-430 signature pad can be used to sign within consult view.

Paired and connected signature pads

None detected.


Pair a Signature Pad

To use a signature pad, you'll need to pair it with your browser.

You should only need to do this once for each device, browser, and signature pad combination.

Pair a signature pad

- d** Then select the STU-430, and click 'Connect'



How to pair signature pad with Concentric

Step 2 - Activate

- e** Next you need to activate the signature pad with a patient. Select a patient and navigate to the “Sign Consent Form” section.
- f** Click on ‘Activate signature pad’

Add your signature to give consent

Sign on this device Sign on a different device **Sign with a signature pad**

Activate your signature pad

A signature pad is plugged in, but you need to activate it before you can use it to sign. Click the button and select the appropriate signature pad from the popup window.

Activate signature pad



- g** You will receive this message and your pad is ready to use for all patients

Add your signature to give consent

Sign on this device Sign on a different device **Sign with a signature pad**

✓ Your signature pad is ready to use

Please use the stylus to add your signature. Tap the ‘I consent’ button on the signature pad when you’re done.



Once you have paired a signature pad with your Concentric account, you will be able to use it with any workstation.