


'Concentric Patient Viewer'

How to open Concentric via ppm+



1

Search for the patient in ppm+.

2

Select the 'Concentric Patient Viewer' on the left-hand side of the ppm+ screen.

3

Concentric log-in screen will appear. Log in to Concentric.

After your first login, the system will remember your details for future access.


The patient's details will be automatically populated into Concentric, without the need to search them.

4

Create the consent episode.

5

Once complete either log out of ppm+, or search for the next patient on ppm+. Concentric will auto close the previous patient session when you re-search.



If accessing Concentric on a kiosk workstation, the browser must be closed, by clicking on the X at the top righthand side of the screen, after each use to prevent inadvertent access under a previous user's session

For further information please contact:

✉ leedsth-tr.digitalconsentsupport@nhs.net